

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-JUNV GP1 glycoprotein Antibody (JUN1)
PRODUCT REFERENCE	ARO-A13164
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Pre-glycoprotein polyprotein GP complex

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Pre-glycoprotein polyprotein GP complex, Pre-GP-C, Stable signal peptide, SSP, Glycoprotein G1, GP1, Glycoprotein G2, GP2, GPC