

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Clostridium botulinum botA/BoNT Antibody (BNTA-4.1)
PRODUCT REFERENCE	ARO-A13167
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Botulinum neurotoxin type A

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Botulinum neurotoxin type A, BoNT/A, Bontoxilysin-A, BOTOX, Botulinum neurotoxin type A1, Botulinum neurotoxin A light chain, LC, 3.4.24.69, Botulinum neurotoxin A heavy chain, HC, botA, atx, bonT