

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human PTP4A1/PRL1 Antibody (Iv0169)
PRODUCT REFERENCE	ARO-A13177
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	PRL-1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-PRL-1, PTP(CAAXI), PTPCAAX1, PTP4A1, Protein tyrosine phosphatase type IVA 1, PRL1, Protein-tyrosine phosphatase of regenerating liver 1, Protein-tyrosine phosphatase 4a1