

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human EGFL7 Antibody (Iv0163)
PRODUCT REFERENCE	ARO-A13200
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	NOTCH4-like protein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-NOTCH4-like protein, Zneu1, Multiple EGF-like domains protein 7, Epidermal growth factor-like protein 7, MEGF7, EGFL7, VE-statin, EGF-like protein 7, Multiple epidermal growth factor-like domains protein 7, Vascular endothelial statin