

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Brazilian scorpion Tityustoxin VII/Toxin VII/Ts7 Antibody (Iv0135)
PRODUCT REFERENCE	ARO-A13233
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Beta-mammal/insect toxin Ts1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Beta-mammal/insect toxin Ts1, Beta-like toxin Ts1, PT-Mice-Ins-beta NaTx6.1, Tityustoxin VII, Toxin VII, Ts VII, Ts7, TsTX-VII, Toxin II-11, Toxin III-10, Toxin T2-IV, Toxin gamma, T-gamma, TiTx gamma, Toxin-g, TsTX-I