

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human CD307e/FCRL5 Antibody (lv0127)
PRODUCT REFERENCE	ARO-A13245
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	FcRL5

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-FcRL5, FCRH5, FCRL5, Fc receptor homolog 5, FcRH5, CD307e, FcR-like protein 5, IRTA2, Immune receptor translocation-associated protein 2, BXMAS1, Fc receptor-like protein 5