



# ProteoGenix

**ProteoGenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

**ProteoGenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

products@proteogenix.fr

## PRODUCT DATASHEET

### OVERVIEW

---

**PRODUCT NAME** InVivoMAb Anti-HIV1 env/Env polyprotein Antibody (Iv0115)

**PRODUCT REFERENCE** ARO-A13260

**CLONALITY** Monoclonal Antibody

**ISOTYPE**

**TARGET** Envelope glycoprotein gp160

### PRODUCT INFORMATION

---

**ORIGIN SPECIES** *Mus musculus (Mouse)*

**FORM** Liquid

**BUFFER** 0.01M PBS (pH 7.4)

**STABILITY & STORAGE** Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

**ALIAS** Anti-Envelope glycoprotein gp160, Env polyprotein, Surface protein gp120, SU, Glycoprotein 120, gp120, Transmembrane protein gp41, TM, Glycoprotein 41, gp41, env, Human Immunodeficiency Virus (HIV), HIV1, HIV-1

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**