

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	InVivoMAb Anti-Human CA9 (Iv0098)
<b>PRODUCT REFERENCE</b>	ARO-A13277
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	Membrane antigen MN

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Membrane antigen MN, P54/58N, Carbonic anhydrase 9, Renal cell carcinoma-associated antigen G250, CA-IX, CAIX, pMW1, Carbonic anhydrase IX, RCC-associated antigen G250, Carbonate dehydratase IX, CA9, MN, G250