

## PRODUCT DATASHEET

### OVERVIEW

---

**PRODUCT NAME** InVivoMAb Anti-Human CXCL12/SDF-1 (lv0090)

**PRODUCT REFERENCE** ARO-A13285

**CLONALITY** Monoclonal Antibody

**ISOTYPE**

**TARGET** Interchrine reduced in hepatomas

### PRODUCT INFORMATION

---

**ORIGIN SPECIES** *Mus musculus (Mouse)*

**FORM** Liquid

**BUFFER** 0.01M PBS (pH 7.4)

**STABILITY & STORAGE** Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

**ALIAS** Anti-Interchrine reduced in hepatomas, hSDF-1, Stromal cell-derived factor 1, SDF-1, Pre-B cell growth-stimulating factor, SDF1, SDF1A, C-X-C motif chemokine 12, SDF1B, CXCL12, IRH, PBSF, hIRH