

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human CCL3/MIP-1-alpha (lv0078)
PRODUCT REFERENCE	ARO-A13296
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Tonsillar lymphocyte LD78 alpha protein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Tonsillar lymphocyte LD78 alpha protein, CCL3, G0S19-1, MIP1A, Macrophage inflammatory protein 1-alpha, LD78-alpha(4-69), SCYA3, SIS-beta, Small-inducible cytokine A3, C-C motif chemokine 3, G0/G1 switch regulatory protein 19-1, PAT 464.1, MIP-1-alph