



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

**PRODUCT NAME** InVivoMAb Anti-Human PDGFD (Iv0074)

**PRODUCT REFERENCE** ARO-A13300

**CLONALITY** Monoclonal Antibody

**ISOTYPE**

**TARGET** IEGF

### PRODUCT INFORMATION

---

**ORIGIN SPECIES** *Mus musculus (Mouse)*

**FORM** Liquid

**BUFFER** 0.01M PBS (pH 7.4)

**STABILITY & STORAGE** Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

**ALIAS** Anti-IEGF, PDGFD latent form, PDGF-D, Spinal cord-derived growth factor B, SCDGF-B, Iris-expressed growth factor, PDGFD, Platelet-derived growth factor D, PDGFD receptor-binding form, SCDGFB

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**