

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human IL18/IL1F4 (Iv0030)
PRODUCT REFERENCE	ARO-A13349
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	IL1F4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-IL1F4, IFN-gamma-inducing factor, Interleukin-1 gamma, IL-18, IL-1 gamma, Interleukin-18, Interferon gamma-inducing factor, IGIF, IL18, Iboctadekin