

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human IL8/CXCL8 (Iv0023)
PRODUCT REFERENCE	ARO-A13356
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	GCP/IL-8 protein II

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-GCP/IL-8 protein II, MDNCF, IL8/NAP1 form III, IL-8(1-77), Neutrophil-activating factor, GCP/IL-8 protein I, GCP/IL-8 protein VI, Interleukin-8, Chemokine (C-X-C motif) ligand 8, Neutrophil-activating protein 1, Monocyte-derived neutrophil chemotacti