

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Human CD120b/TNFRSF1B/TNFR2 (SAA0309)
PRODUCT REFERENCE	ARO-A13374
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	TNFR2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-TNFR2, TNFR-II, p80 TNF-alpha receptor, Etanercept, TNFBR, TBP-2, TBPII, Tumor necrosis factor receptor superfamily member 1B, CD120b, p75, TNFRSF1B, Tumor necrosis factor receptor 2, TNF-RII, TNF-R2, Tumor necrosis factor receptor type II