

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-Human LCP1 Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A13424 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | E. coli - derived recombinant Human LCP1 (Thr106-Lys379). |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-Lymphocyte cytosolic protein 1, LC64P, L-plastin, LCP-1, Plastin-2, LCP1, PLS2, |