

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human EIF3H Polyclonal Antibody
PRODUCT REFERENCE	ARO-A13437
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human EIF3H (Lys37-Ser208).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-EIF3H, eIF-3-gamma, eIF3h, eIF3 p40 subunit, EIF3S3, Eukaryotic translation initiation factor 3 subunit 3, Eukaryotic translation initiation factor 3 subunit H