

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-Human KLRG1/CLEC15A Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A13455 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | E. coli - derived recombinant Human KLRG1/CLEC15A (Ser74-Lys186). |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-Killer cell lectin-like receptor subfamily G member 1, C-type lectin domain family 15 member A, ITIM-containing receptor MAFA-L, MAFA-like receptor, Mast cell function-associated antigen, KLRG1, CLEC15A, MAFA, MAFAL |