

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Coat protein Polyclonal antibody
PRODUCT REFERENCE	ARO-A13508
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant ChiVMV Coat protein (Ser2816-Leu3089).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-49 kDa proteinase,6 kDa protein 1,6 kDa protein 2,Coat protein,Cytoplasmic inclusion protein,Genome polyprotein,Helper component proteinase,N-terminal protein,N1a-pro,Nuclear inclusion protein A,Nuclear inclusion protein B,P1 proteinase,RNA-directed