

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ELOVL3/CIG30 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13509
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Mouse ELOVL3 / CIG30 (Asp110-Gly169).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-3-keto acyl-CoA synthase ELOVL3,ELOVL FA elongase 3,Elongation of very long chain fatty acids protein 3,CIG30,ELOVL3,Cold-inducible glycoprotein of 30 kDa,Very long chain 3-oxoacyl-CoA synthase 3,ELOVL fatty acid elongase 3,Very long chain 3-ketoacyl