

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SCAF11 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13516
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human SCAF11 (Arg944-Trp1148).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-SFRS2IP, SR-related and CTD-associated factor 11, SIP1, SRSF2IP, Splicing regulatory protein 129, CASP11, Splicing factor, arginine/serine-rich 2-interacting protein, CTD-associated SR protein 11, Renal carcinoma antigen NY-REN-40, SC35-interacting protein 1, S