

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Adiponectin/ACRP30/ADIPOQ Polyclonal antibody
<b>PRODUCT REFERENCE</b>	ARO-A13519
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	E. coli - derived recombinant Human Adiponectin / ACRP30 / ADIPOQ (Glu19-Asn244).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Adiponectin,ACRP30,GBP28,Adipose most abundant gene transcript 1 protein,ACDC,APM1,ADIPOQ,apM-1,30 kDa adipocyte complement-related protein,Adipocyte complement-related 30 kDa protein,Adipocyte,C1q and collagen domain-containing protein,Gelatin-bind