

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-FUT2 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13522
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human FUT2 (Gln29-His343).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Fucosyltransferase 2, SE2, Se, Galactoside alpha-(1,2)-fucosyltransferase 2, Secretor blood group alpha-2-fucosyltransferase, Alpha(1,2)FT 2, Type 1 galactoside alpha-(1,2)-fucosyltransferase FUT2, FUT2, GDP-L-fucose:beta-D-galactoside 2-alpha-L-fucosyltrans