

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-S100A12 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13523
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human S100A12 (Met1-Glu92).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Migration inhibitory factor-related protein 6,MRP-6,Extracellular newly identified RAGE-binding protein,p6,Calcium-binding protein in amniotic fluid 1,Calgranulin-C,Protein S100-A12,CGRP,S100 calcium-binding protein A12,Neutrophil S100 protein,CAAF1,