

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-S6K1/RPS6KB1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13530
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human S6K1 / RPS6KB1 (Phe91-Phe352).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-S6K-beta-1,70 kDa ribosomal protein S6 kinase 1,p70 S6KA,p70-S6K 1,S6K1,Serine/threonine-protein kinase 14A,p70 ribosomal S6 kinase alpha,p70 S6 kinase alpha,STK14A,RPS6KB1,Ribosomal protein S6 kinase I,p70 S6K- alpha,P70S6K1,Ribosomal protein S6 kina