

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PTN Polyclonal antibody
PRODUCT REFERENCE	ARO-A13531
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human PTN (Gly33-Asp168).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Heparin-binding neurite outgrowth-promoting factor 1, Heparin-binding brain mitogen, PTN, HBGF-8, Pleiotrophin, Heparin-binding growth factor 8, HBBM, HBNF-1, NEGF1, Osteoblast-specific factor 1, HBNF, HBNF1, Heparin-binding growth-associated molecule, HB-GAM, OSF