

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-XRCC5/Ku80 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13533
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human XRCC5 / Ku80 (Leu384-Ile732).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Lupus Ku autoantigen protein p86,XRCC5,Thyroid-lupus autoantigen,Ku86,ATP-dependent DNA helicase 2 subunit 2,CTC85,G22P2,DNA repair protein XRCC5,86 kDa subunit of Ku antigen,TLAA,Nuclear factor IV,CTC box-binding factor 85 kDa subunit,ATP-dependent