

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PPBP/CXCL7/NAP-2 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13535
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human PPBP / CXCL7 / NAP-2 (Ala59-Asp128).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-NAP-2,NAP-2(1-63),NAP-2(73),PPBP,CXCL7,PBP,Small-inducible cytokine B7,Macrophage-derived growth factor,LA-PF4,Leukocyte-derived growth factor,CTAP-III,CTAP-III(1-81),Low-affinity platelet factor IV,TGB1,Beta-TG,C-X-C motif chemokine 7,MDGF,Platelet