

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ECI2 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13540
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human ECI2 (Met31-Val241).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Renal carcinoma antigen NY-REN-1,PECI,DRS-1,ECI2,Dodecenoyl-CoA isomerase,D3,D2-enoyl-CoA isomerase,HCA88,Peroxisomal 3,2-trans-enoyl-CoA isomerase,Hepatocellular carcinoma-associated antigen 88,Enoyl-CoA delta isomerase 2,pECI,DRS1,DBI-related prote