

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ARNTL Polyclonal antibody
PRODUCT REFERENCE	ARO-A13541
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human ARNTL (Glu48-Leu320).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Basic-helix-loop-helix-PAS protein MOP3,ARNTL,Member of PAS protein 3,PAS domain-containing protein 3,BMAL1,Aryl hydrocarbon receptor nuclear translocator-like protein 1,BHLHE5,Class E basic helix-loop-helix protein 5,bHLHe5,MOP3,Brain and muscle ARN