

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-FKBP7 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13568
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human FKBP7 (Gln24-Tyr217).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-23 kDa FK506-binding protein,FKBP-7,FKBP23,FKBP7,Rotamase,PPlase FKBP7,FK506-binding protein 7,Peptidyl-prolyl cis-trans isomerase FKBP7,FKBP-23,23 kDa FKBP,