



# ProteoGenix

## ProteoGenix S.A.S

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## ProteoGenix INC.

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

|                   |                                                          |
|-------------------|----------------------------------------------------------|
| PRODUCT NAME      | Anti-PSAT1 Polyclonal antibody                           |
| PRODUCT REFERENCE | ARO-A13571                                               |
| CLONALITY         | Polyclonal Antibody                                      |
| ISOTYPE           | IgG                                                      |
| TARGET            | E. coli - derived recombinant Human PSAT1 (Met1-Gly312). |

### PRODUCT INFORMATION

---

|                     |                                                                                                                                                                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES      | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                               |
| FORM                | Liquid                                                                                                                                                                                                                                                                                                              |
| BUFFER              | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300                                                                                                                                                                                                                                                                 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS               | Anti-Phosphohydroxythreonine aminotransferase,PSAT,PSAT1,PSA,Phosphoserine aminotransferase,                                                                                                                                                                                                                        |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY