

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-DIDO1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13665
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human DIDO1 (Lys258-Cys455).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Death-inducer obliterator 1,DIO- 1,DATF1,C20orf158,KIAA0333,hDido1,DIDO1,Death-associated transcription factor 1,DATF-1,