

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-IL31RA Polyclonal antibody
PRODUCT REFERENCE	ARO-A13715
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human IL31RA (Val130-Asp622).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-IL-31RA, Interleukin-31 receptor subunit alpha, IL31RA, hGLM-R, ZcytoR17, IL-31R-alpha, GPL, CRL3, Cytokine receptor-like 3, IL-31R subunit alpha, IL-31 receptor subunit alpha, Gp130-like receptor, GLM-R, Gp130-like monocyte receptor,