

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-SND1 Polyclonal antibody |
| PRODUCT REFERENCE | ARO-A13742 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | E. coli - derived recombinant Human SND1 (Ala650-Arg910). |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-TDRD11, SND1, p100 co-activator, Staphylococcal nuclease domain-containing protein 1, Tudor domain-containing protein 11, EBNA2 coactivator p100, 100 kDa coactivator, |