

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-H3C15 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13744
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human H3C15 (Met1-Ala136).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HIST2H3A,H3C15,Histone H3/o,H3-clustered histone 14,H3-clustered histone 15,H3-clustered histone 13,H3F2,H3FM,Histone H3/m,HIST2H3D,HIST2H3C,Histone H3.2,