

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SMAD1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13786
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human SMAD1 (Asp259-Ser462).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-BSP1, Mothers against DPP homolog 1, Mothers against decapentaplegic homolog 1, SMAD1, Transforming growth factor-beta-signaling protein 1, MADH1, JV4-1, MADR1, SMAD family member 1, hSMAD1, SMAD 1, MAD homolog 1, Mad-related protein 1, BSP-1, Smad1,