

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-MAPK11 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13787
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human MAPK11 (Leu24-Phe308).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-MAP kinase 11, Mitogen-activated protein kinase p38 beta, MAPK 11, Mitogen-activated protein kinase 11, PRKM11, SAPK2B, Stress-activated protein kinase 2b, MAP kinase p38 beta, MAPK11, p38b, p38-2, SAPK2, SAPK2b,