

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |                                                             |
|--------------------------|-------------------------------------------------------------|
| <b>PRODUCT NAME</b>      | Anti-MMRN1 Polyclonal antibody                              |
| <b>PRODUCT REFERENCE</b> | ARO-A13829                                                  |
| <b>CLONALITY</b>         | Polyclonal Antibody                                         |
| <b>ISOTYPE</b>           | IgG                                                         |
| <b>TARGET</b>            | E. coli - derived recombinant Human MMRN1 (Ala807-Gly1053). |

### PRODUCT INFORMATION

---

|                                |                                                                                                                                                                                                                                                                                                                     |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORIGIN SPECIES</b>          | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                               |
| <b>FORM</b>                    | Liquid                                                                                                                                                                                                                                                                                                              |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300                                                                                                                                                                                                                                                                 |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | Anti-MMRN,ECM,Multimerin-1,EMILIN4,p155,EMILIN-4,p-155,Endothelial cell multimerin,GPIA*,MMRN1,Elastin microfibril interface located protein 4,Elastin microfibril interfacier 4,                                                                                                                                   |