

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PRKAA1/AMPK alpha 1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13832
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human PRKAA1 / AMPK alpha 1 (Met10-Gln559).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-AMPK subunit alpha-1, Acetyl-CoA carboxylase kinase, Tau-protein kinase PRKAA1, ACACA kinase, Hydroxymethylglutaryl-CoA reductase kinase, 5'-AMP-activated protein kinase catalytic subunit alpha-1, AMPK1, PRKAA1, HMGCR kinase,