

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-POU5F1/OCT3/OCT4 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13857
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human POU5F1 / OCT3 / OCT4 (Ser136-Ser289).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Octamer-binding protein 4,Oct-4,OCT3,OCT4,OTF3,OTF-3,Octamer-binding transcription factor 3,POU5F1,Oct-3,POU domain,class 5,transcription factor 1,Octamer-binding protein 3,