

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-OLR1 Polyclonal antibody
<b>PRODUCT REFERENCE</b>	ARO-A13872
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	E. coli - derived recombinant Human OLR1 (Leu81-Gln273).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-C-type lectin domain family 8 member A, LOX-1, Ox-LDL receptor 1, Lectin-like oxidized LDL receptor 1, OLR1, Oxidized low-density lipoprotein receptor 1, Lectin-like oxLDL receptor 1, LOX1, CLEC8A, hLOX-1, Lectin-type oxidized LDL receptor 1,