

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SEC11A Polyclonal antibody
PRODUCT REFERENCE	ARO-A13878
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human SEC11A (Met1-Glu179).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Microsomal signal peptidase 18 kDa subunit,Endopeptidase SP18,SPCS4A,Signal peptidase complex catalytic subunit SEC11A,SEC11A,SPase 18 kDa subunit,SEC11-like protein 1,SEC11 homolog A,SPC18,SEC11L1,