

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-UBE2I/UBC9 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13880
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human UBE2I / UBC9 (Asp33-Ser158).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-UBE2I,UBC9,p18,UBCE9,Ubiquitin carrier protein 9,RING-type E3 SUMO transferase UBC9,SUMO-conjugating enzyme UBC9,Ubiquitin carrier protein I,Ubiquitin-protein ligase I,Ubiquitin-conjugating enzyme E2 I,SUMO-protein ligase,