

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SUMO1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13882
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human SUMO1 (Met1-Val101).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-SUMO-1, SMT3C, SMT3H3, Sentrin, GMP1, Ubiquitin-like protein UBL1, UBL1, Small ubiquitin-related modifier 1, Ubiquitin-homology domain protein PIC1, SUMO1, Ubiquitin-like protein SMT3C, SMT3 homolog 3, Smt3C, GAP-modifying protein 1,