

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CCL23/MIP-3/MPIF-1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13908
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human CCL23 / MIP-3 / MPIF-1 (Arg22-Asn120).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-CKB-8,MIP3,CCL23,C-C motif chemokine 23,MPIF1,MIP-3,Macrophage inflammatory protein 3,SCYA23,Small-inducible cytokine A23,Myeloid progenitor inhibitory factor 1,CK-beta-8,MPIF-1,