

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-MRE11 Polyclonal antibody
PRODUCT REFERENCE	ARO-A13942
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human MRE11 (Pro375-Lys609).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HNGS1, Double-strand break repair protein MRE11, MRE11A, Double-strand break repair protein MRE11A, Meiotic recombination 11 homolog A, MRE11 homolog A, MRE11 homolog 1, MRE11, Meiotic recombination 11 homolog 1,