

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-KRT6C Polyclonal antibody
PRODUCT REFERENCE	ARO-A13949
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human KRT6C (Glu163-Gly482).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Type-II keratin Kb12,KRT6E,CK-6C,Keratin K6h,Keratin-6C,Cytokeratin-6E,Keratin,type II cytoskeletal 6C,K6C,KRT6C,Cytokeratin-6C,CK-6E,