

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |                                                             |
|--------------------------|-------------------------------------------------------------|
| <b>PRODUCT NAME</b>      | Anti-SERPINB6 Polyclonal antibody                           |
| <b>PRODUCT REFERENCE</b> | ARO-A13999                                                  |
| <b>CLONALITY</b>         | Polyclonal Antibody                                         |
| <b>ISOTYPE</b>           | IgG                                                         |
| <b>TARGET</b>            | E. coli - derived recombinant Human SERPINB6 (Met1-Pro376). |

### PRODUCT INFORMATION

---

|                                |                                                                                                                                                                                                                                                                                                                     |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORIGIN SPECIES</b>          | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                               |
| <b>FORM</b>                    | Liquid                                                                                                                                                                                                                                                                                                              |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300                                                                                                                                                                                                                                                                 |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | Anti-CAP,PI6,PTI,PI-6,Cytoplasmic antiproteinase,Serpin B6,Peptidase inhibitor<br>6,SERPINB6,Placental thrombin inhibitor,                                                                                                                                                                                          |