

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-GDI1 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14010
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human GDI1 (Met1-Val274).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Oligophrenin-2,Rab GDP dissociation inhibitor alpha,XAP4,Guanosine diphosphate dissociation inhibitor 1,GDIL,GDI1,RABGDIA,Rab GDI alpha,Protein XAP-4,OPHN2,GDI-1,