

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PTPN6 Polyclonal antibody
PRODUCT REFERENCE	ARO-A14023
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human PTPN6 (Phe244-Phe515).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HCP,PTPN6,PTP1C,Hematopoietic cell protein-tyrosine phosphatase,SH-PTP1,Protein-tyrosine phosphatase SHP-1,Protein-tyrosine phosphatase 1C,Tyrosine-protein phosphatase non-receptor type 6,PTP-1C,